

Galaxy Tutorial - Basics of Data Analysis using Galaxy platform, Condensed

(SW4) The Galaxy Platform for Multi-Omic Data Analysis and Informatics

Dave Clements, Galaxy Project, Johns Hopkins University
8:30-10:00am, Saturday February 20, 2016

ABRF 2016

[Explore Galaxy](#)

[Left/Tools panel](#)

[Right/History panel](#)

[Center panel](#)

[Create an account](#)

[Our goal: Identify transcripts using RNA-Seq in pre-pro-B and pro-B mouse cells](#)

[The data](#)

[Get the data / Data libraries](#)

[FASTQ datasets](#)

[FASTQ Format](#)

[Check FASTQ dataset quality](#)

[Run our first tool: FastQC](#)

[Per base sequence quality \(boxplot\)](#)

[Per tile sequence quality](#)

[Per sequence quality scores](#)

[Per base sequence content](#)

[Adapter Content](#)

[And the rest](#)

[Rerun FastQC on pro-B dataset](#)

[Improve the quality: Trimmomatic](#)

[Processing Multiple datasets with one submission](#)

[Select Trimmomatic Filters](#)

[Trimmomatic Results](#)

[Using Galaxy Scratchpad to view multiple datasets](#)

[Map the reads](#)

[HISAT2](#)

[HISAT2 Results](#)

[Updating dataset and history metadata: Renaming](#)

[But what do the HISAT2 mappings look like?](#)

[Transcript Prediction with StringTie](#)

[And there is more!](#)

[Workflows: repeating an analysis](#)

[Working with multiple histories](#)


[Sharing and publishing](#)

[Galaxy pages](#)

A longer version of the same material is available online at <http://bit.ly/gxyabrf16intro-long>

This handout is available online at <http://bit.ly/gxyabrf16intro>

Create an account



To create an account, *click* on the **User** pulldown and **select Register**.

Use an email address you can remember, **and a low security password**. Note that these servers use HTTP, not HTTPS, and therefore your password is going out on the net unencrypted. (Galaxy can use HTTPS and most servers do, but the default cloud installation uses HTTP.)

A screenshot of a web form titled 'Create account'. It contains the following fields: 'Email address' (with the value 'clements@galaxyproject.org'), 'Password' (with masked characters '*****'), 'Confirm password' (with masked characters '*****'), and 'Public name' (with the value 'clements'). Below the fields is a 'Submit' button. A note at the bottom states: 'Your public name is an identifier that will be used to generate addresses for information you share publicly. Public names must be at least three characters in length and contain only lower-case letters, numbers, and the '-' character.'

Reload the page by *clicking* **Analyze Data** or **Return to the home page**.

Our goal: Identify transcripts using RNA-Seq in pre-pro-B and pro-B mouse cells

We are going to use RNA-seq data to detect transcripts, possibly novel transcripts, in pre-pro-B and pro-B mouse cells.


This data comes from

Heydarian *et al.*, (2014) Prediction of Gene Activity in Early B Cell Development Based on an Integrative Multi-Omics Analysis. *J Proteomics Bioinform* 7: 050-063. doi:10.4172/jpb.1000302

The paper uses RNA-Seq, ChIP-Seq, GRO-Seq, and iTRAQ techniques to “demonstrate that active chromatin modifications at promoters are good indicators of transcription and steady state mRNA levels.”

The data

We have two RNA-Seq datasets, one from mouse pre-pro-B cells, and one from mouse pro-B cells. Each dataset consists of ~90 million single-end reads. Sequencing was done on an Illumina HiSeq 2000 and all reads are 97bp long.



We are going to use stripped down datasets consisting of reads that tend to map to the particular regions of the genome we'll focus on today, plus some additional reads, just to make things interesting.



Note that while reducing the datasets will make tool execution very fast, it can be a dangerous cheat. Many NGS tools are built on statistical models that assume they are being given data for the entire genome. With our data that is decidedly not the case. Downsampling your data is still a useful technique for early experimentation when building your analysis, but results and assumptions should be checked when scaling up your analysis to complete datasets.

Get the data / Data libraries

Our reduced datasets are available in a *data library* on your server. Click on the **Shared Data** pull-down, and then select **Data Libraries**.

Then click on the **RNA-Seq Transcript Prediction** library and select all 3 datasets in this library and then click on **to History** to import them into your current history.


This will ask which history to import these datasets into. So far, you only have one ("Unnamed history"). Click **Import** to add them to that history.



Now click **Analyze Data** to see what we have.


FASTQ datasets

And what we have is a history consisting of three datasets. Two of them are the raw read datasets for the pre-pro-B and pro-B cells. These datasets are in FASTQ format.



To see a preview of any dataset, click on the dataset name. This displays some metadata about the dataset, and if it's in a text-based format, it will also show the first few lines of the dataset.


To view the full dataset, *poke it in the eye*. This will display the content of this FASTQ dataset in the middle panel.



FASTQ Format

Each FASTQ entry is 4 lines long and incorporates three types of information:

- Line 1 starts with an @ and is the read's identifier
- Line 2 shows the called bases in this read
- Line 3 is a separator (the + character).
- Line 4 shows the instrument's confidence in each of the base calls.



https://en.wikipedia.org/wiki/FASTQ_format

Phred quality scores are assigned to each base call. Phred scores are logarithmic and represent the likelihood that given base call is incorrect. Most datasets use Sanger scoring which typically ranges from 0 to 40. Some reference points in that spectrum:

Score	Instrument thinks it will be wrong	
10	10% of the time	(1 out of 10 times)
20	1% of the time	(1 out of 100 times)
30	0.1% of the time	(1 out of 1000 times)
40	0.01% of the time	(1 out of 10000 times)

We don't see these scores in the FASTQ files because each score is encoded to a single character according the FASTQSanger convention. FASTQSanger maps scores to a range of adjacent ASCII characters, starting with ! for 0, and ending with I for 40.

Check FASTQ dataset quality

There are two ways to find tools on a Galaxy server. If you know the name of the tool you can enter it in the tool panel search box. You can also try searching for more general terms like “trim” or “quality.” Searching with more general terms will be hit or miss.

You can also just browse each Toolbox. In this case we are investigating the quality of an NGS dataset, so the **NGS: QC and Manipulation** toolbox looks the most promising. *Click* on the toolbox name to see the tools in it. In this case **FastQC Read Quality reports** is right at the top, and that sounds exactly like what we are looking for. *Click* on it.

FastQC is a widely used tool for summarizing the quality of FASTQ datasets. It's good at presenting the big picture, and at identifying specific

[Metagenomic analyses](#)
[Motif Tools](#)
[NGS: QC and manipulation](#)
[FastQC Read Quality reports](#)
[Tabular to FASTQ converter](#)
[FASTQ Quality Trimmer by sliding window](#)
[FASTQ Trimmer by column](#)

areas that may need to be addressed.

This is a typical, if simple, Galaxy tool form. Tool pages generally have these sections

- *Versions and options* - If multiple versions are supported, there will be link to run the others. And the options menu enables you see more about the tool and how it is wrapped in Galaxy.
- *The form* - where you configure and run the tool
- *Documentation* - Describes the tool and provides help. This section often include links to a tool's external documentation.
- *Citation* - A paper, if there is one, for this tool.












Run our first tool: FastQC

FastQC has only three parameters that can be set. Let's set only one, **Short read data from your current history**. *Click* this option's dataset pulldown menu and select **2: pre-pro-B chrX/12 regions FASTQ raw**. We aren't going to set the contaminant list or limit the submodules/tests that FastQC runs. We'll request that it run all of its tests. *Click* **Execute**.

Now a quick succession of things (hopefully) happens

1. A **big green box** appears in the middle panel and two **small gray boxes** appear at the top of the history. These mean that the task has been successfully queued: It's ready to run, but isn't running yet.
2. The gray boxes in the history are replaced with **yellow boxes**. This means the job has started and is actively running
3. The yellow box are replaced with **green boxes** indicating that the tool finished successfully.

Summary

-  [Basic Statistics](#)
-  [Per base sequence quality](#)
-  [Per tile sequence quality](#)
-  [Per sequence quality scores](#)
-  [Per base sequence content](#)
-  [Per sequence GC content](#)
-  [Per base N content](#)
-  [Sequence Length Distribution](#)
-  [Sequence Duplication Levels](#)
-  [Overrepresented sequences](#)
-  [Adapter Content](#)
-  [Kmer Content](#)

Once a tool has finished running you can preview the output datasets (by clicking on the dataset name) or view the data itself (by poking it in the eye). **Poke the FastQC Web page dataset in the eye.**


This brings up 12 different reports about the pre-pro-B FASTQ dataset.

For each test it runs, FastQC has predefined thresholds for what constitutes good, not-so-good, and bad. This dataset has passed 7 of the tests and got warning on the rest. Let's take a look at some of these tests.

Per base sequence quality (boxplot)

Boxplots are the most common way to summarize the quality of a FASTQ dataset.

- **X axis is position in the read.** That is, the first column summarizes all the base calls in the first position of all the reads; the second column summarizes the second call in all the reads and so on. Starting at position 10, the information is aggregated for every two bases.
- **Y axis is Phred quality score**, from 0 to 40.
- Each column summarizes the quality scores across all reads at that position
 - **Yellow boxes** bound the middle 50% of quality scores at each position in the read.
 - **Whiskers** bound the middle 80% of quality scores at that position.
 - **Blue line** is the average score at each position
 - **Red lines** are the median score at each position.





Per tile sequence quality

This graph is particularly significant for this crowd. It highlights where on the slide the low quality reads tended to be. As core staff you can use this information to identify trouble spots with the sequencing run.

For Illumina data, tile information is embedded in the ID line of FASTQ entries.

Note: We're showing the graphic for the complete ~90 million read dataset.





Per sequence quality scores

This graph gives us an idea just how bad (or good) our bottom 10% is. Here, the X axis is the quality score, and the Y axis is the number of times that score occurs across the entire dataset. The scores in this dataset are pretty good and there are only a small number of calls with confidence below 20. However, there is a small spike at 2. We'll want to drop those.



Per base sequence content

This graph is good at highlighting untrimmed adapter or barcode problems.

Flat lines mean there is no bias for particular bases at particular positions. The turbulence at the front is caused by the use of random hexamers in Illumina library preparation. See

Hansen, *et al.*, "[Biases in Illumina transcriptome sequencing caused by random hexamer priming](#)" *Nucleic Acids Research*, Volume 38, Issue 12 (2010)

The summary is: *It's not our fault.*



Adapter Content

At the end of the reads the dataset is 3 to 4 percent untrimmed adapter. This falls within FastQC's tolerance, but not ours.

And the rest


The remaining tests are either not informative (look at Sequence Length Distribution), flat out good (Per base N content), or look significantly different when the whole dataset is scanned (Per sequence GC content, Sequence Duplication Levels, Kmer Content). We'll ignore all these today.

Rerun FastQC on pro-B dataset

Now let's run FastQC on the other dataset, the pro-B reads. We could rerun FastQC by finding it in the Tool panel again, or, we could rerun FastQC from the History panel. Let's do that.

Preview one of the **FastQC** datasets in your history by *clicking* on its name. In the preview *click* on the **Run this job again** looping arrow icon.

This launches the tool form in the middle panel and pre-populates it with the same settings it was run with before.

A screenshot of the FastQC tool form in Galaxy. The title bar says "FastQC Read Quality reports (Galaxy Tool Version 0.63)". Under "Short read data from your current history", a dropdown menu shows "2: pre-pro-B chrX/12 regions FASTQ raw" selected. Below, the "Contaminant list" is empty. Another dropdown shows "3: pro-B chrX/12 regions FASTQ raw" selected. The "Submodule and Limit specifying file" dropdown is set to "Nothing selected". An "Execute" button is at the bottom left.

Click the dataset pulldown menu and *select* **3: pro-B chrX/12 regions FASTQ raw**. Then *click* **Execute**. *Poke* the new **FastQC Webpage** in the eye to review the results. The quality is similar.



Improve the quality: Trimmomatic

Use Trimmomatic to separate the good from the bad.

FastQC was the first tool in the **NGS: QC and manipulation** toolbox. **Trimmomatic** is the last. Click on it. First, *change* **Paired end data?** to **No**.

Processing Multiple datasets with one submission

We could run Trimmomatic the same way we ran FastQC - by running it on one dataset, and then re-running it on the 2nd.

A screenshot of the Trimmomatic tool form in Galaxy. The title bar says "Trimmomatic flexible read trimming tool for Illumina NGS data (Galaxy Tool Version 0.32.3)". Under "Paired end data?", the "No" button is selected. The "Input FASTQ file" dropdown shows "3: pro-B chrX/12 regions FASTQ raw" selected.A screenshot of the Trimmomatic tool form in Galaxy, similar to the previous one. A mouse cursor is clicking on the "Multiple datasets" button, which is located below the "Input FASTQ file" dropdown. The "Paired end data?" section shows "No" selected.

We can achieve the same effect more efficiently by using the **Multiple datasets** feature, another available shortcut.

To enable this *click* on the **Multiple datasets** icon under **Input FASTQ file**.

This changes the parameter selection from a pull-down to a multi-select field. *Select* both datasets.

Select Trimmomatic Filters


We want to use three different filters in Trimmomatic:

1. **ILLUMINACLIP** - This will trim known Illumina adapters from reads in the dataset.
2. **SLIDINGWINDOW** - Cutting once the average quality within a sliding window falls below a threshold. Trims from both ends
3. **MINLEN** - Drop any reads that are shorter than a given length after trimming.


First, *select Yes* under **Perform initial ILLUMINACLIP step?**. Then *change* the **Adapter sequences to use** pulldown to **TruSeq3 (single-ended, for MiSeq and HiSeq)**. Leave other parameters at the default.

Update the default **SLIDINGWINDOW** step to use the same window definition as the paper. *Change Number of bases to average across* to **3**.

To deal with reads that are too short add another operation.



Click on **+ Insert Trimmomatic Operation**. *Change Select Trimmomatic operation to perform* to **Drop reads below a specified length (MINLEN)** and *set Minimum length of reads to be kept* to **32**. All 3 steps are set up. *Click Execute*.




Trimmomatic Results

To see the net decrease in dataset size, preview the before and after datasets, and compare their size.

To see the full Trimmomatic summary, *click* the ⓘ icon (**View details**) in the dataset preview. This displays the metadata for this dataset.

Two links for **stdout** and **stderr** are included with the metadata. These are standard files in Unix to record *standard output* and *standard error* from the tool.


Click stdout to see the full summary. This is informative, but let's run FastQC again on both and compare those reports with the original reports.





Tool version:	
Tool Standard Output:	stdout
Tool Standard Error:	stderr
Tool Exit Code:	0


Click **FastQC** in the tool panel, select the **Multiple datasets** icon under **Short read data** from your current history, and then select both **Trimmomatic on ...** datasets. Click **Execute**.

Note: We aren't using paired-end data today, but if we were Trimmomatic would automatically maintain read pairings throughout all steps.



Using Galaxy Scratchpad to view multiple datasets


We want to compare pre and post-Trimmomatic FastQC reports. Scratchpad will enable us to do this without going insane. To enable Scratchpad, click on the grid icon () in the top menu bar. This turns it yellow ()! Now, scroll down in your history to first FastQC report on the pre-pro-B data (most likely named **4: FastQC on data 2: Webpage**) and poke it in the eye.



Click somewhere outside the displayed dataset. Scroll up the history to the second FastQC report on the pre-pro-B data (most likely the fourth dataset from the top) and poke it in the eye.

After some scrolling and resizing you can now compare the two reports at once.. Additional datasets can be added as desired.

We see that the quality has improved.

To switch back to the standard view either click outside the displayed datasets or poke the yellow eye (). To disable the Scratchbook display, click on the yellow grid icon.

Map the reads

The data has now been cleaned and we are ready to map them against the mouse reference genome.

HISAT2

Bring up the HISAT2 tool form and **change Single end or paired reads? to Individual unpaired reads** and then **select the Multiple datasets icon** under **Reads**. **Select both Trimmomatic datasets**.

HISAT is an “option-rich” tool. The HISAT tool wrapper takes a common approach of placing them in sections that can be expanded as needed. HISAT has 5 such sections:

If you expand all sections, HISAT options go from 1 screen to almost 6. Understanding all these options is the right thing to do, but it’s also daunting. One of Galaxy’s strengths is that it allows you to *experiment with tools* and *learn them incrementally*. We are going to start with setting just one parameter.

Change Spliced alignment parameters to Specify spliced alignment parameters. This causes almost 15 additional options to be displayed, most of them giving you fine-grained control over how reads that are split across multiple exons are scored. Ignore all those and *scroll* down to **GTF file with known splice sites** and **select 1: GTF Annotation chrX/12 regions** and **click Execute**.

HISAT2 Results

HISAT2 creates a mapped reads dataset for each of the two inputs.

These files are in BAM format, an efficient binary format for representing mapped data.

To see the full HISAT2 summary *click* the ① icon (**View details**) in the dataset preview. Then *click* the **stderr** link. Our overall alignment rates are over 98% for both datasets.

HISAT2 A fast and sensitive alignment program (Galaxy Tool Version 1.0.1) Options

Input data format
FASTQ

Single end or paired reads?
Individual unpaired reads

Reads
[Batch mode input field containing 4 datasets]
This is a batch mode input field. A separate job will be triggered for each dataset.

Source for the reference genome to align against
Use a built-in genome
Built-in references were created using default options

Select a reference genome
Mouse Dec. 2011 (GRCm38/mm10) (mm10)
If your genome of interest is not listed, contact the Galaxy team

Alignment options
Use default values

Input options
Use default values

Scoring options
Use default values

Spliced alignment parameters
Use default values

Paired alignment parameters
Use default values


Disable spliced alignment
Yes No

GTF file with known splice sites
[Batch mode input field containing 1 dataset]
1: GTF Annotation chrX/12 regions

Transcriptome assembly reporting
☒ Use default reporting.
☐ Report only those alignments within known transcripts.
☐ Report alignments tailored for transcript assemblers including StringTie.
☐ Report alignments tailored specifically for Cufflinks.

```
[bam_header_read] EOF marker is absent. The input is probably truncated.
[samopen] SAM header is present: 66 sequences.
31843 reads; of these:
  31843 (100.00%) were unpaired; of these:
    482 (1.51%) aligned 0 times
    26747 (84.00%) aligned exactly 1 time
    4614 (14.49%) aligned >1 times
98.49% overall alignment rate
```


Updating dataset and history metadata: Renaming



The HISAT2 output datasets have names like **HISAT2 on data 1 and data 8**. Rename them to clearly indicate what's in them.

Click on the **pencil (Edit attributes)** icon for the pre-pro-B

HISAT2 dataset. Name it something like **HISAT2 pre-pro-B mapped reads**. Click **Save** and then *repeat* with the **pro-B dataset**.



We should also



give our history a name. A best practice is to always name your histories. To do this, *click* on **Unnamed history**, *enter* a new and informative history name, and then *press return* or **enter**.

But what do the HISAT2 mappings look like?


There are summary tools that provide more statistics on things like mapping state (**Flagstat**) and how many reads mapped to each chromosome (**IdxStats**).

We can visually inspect the mappings in the context of the genomic regions they mapped to. Later today you'll use the IGV desktop program to visualize genomic data. Right now, let's learn how to do this inside Galaxy, using the Trackster visualization tool.

To launch Trackster, open a preview of one of the HISAT2 datasets and *click* on the **Visualize in Trackster** icon.





Select **View in new visualization** and then give the visualization a meaningful name and *click* **Create**.




This will bring up a new Trackster browser with a status message saying that it is indexing the dataset for display. *Switch* from chr1 to **chr12** (or chrX), where we have mapped data.




While the track is loading, add other datasets to the display. *Click* on the **+** (**Add tracks**) icon in the upper right, and then select the other mapped dataset and the reference annotation.





Once those new tracks have loaded, zoom in on the region of interest - the area with the annotations: *Click and drag* on the coordinates track (near the top) and include the entire region with annotation.



You may have to repeat that a few times to get the region of interest:




The top 2 tracks show the read depth for any position in this region, and the bottom track shows the annotation (but not very well). *Hover* over the **GTF Annotation chrX/12 regions** track and *click* the **Set display mode icon** and *select* **Squish**. To get a better idea of what the aligned reads look like, zoom in the last 5 or so introns on the far right.



This changes the display mode for all tracks and now we can see individual reads, including where splice junction were introduced to map the reads. Drag the annotation track to the middle by *clicking and dragging* the **vertical lines icon** at the left of the track title upwards to between the two BAM tracks.

We'll want to come back to this location once we have predicted transcripts for these datasets. Bookmark this location by *clicking* on the **bookmark ribbon icon** in the upper right corner. Then *click* the **+ (Add bookmark) icon** and give the bookmark a meaningful description. Finally, *click* the **save icon** to save your visualization.




Transcript Prediction with StringTie

HISAT2 makes no prediction about how those reads might assemble into transcripts. That job is handled by StringTie.

Search for **StringTie** in the tool panel and launch the tool form. Select **Multiple datasets** under **Mapped reads to assemble transcripts from**, and then select both **HISAT2** datasets. Change **Use GFF file to guide assembly** to **Use GFF**. Click **Execute**.


StringTie generates 3 output datasets for each input dataset:

- Coverage - All the transcripts in the reference GTF that are fully covered by reads.
- Gene abundance estimates - This is created but it's empty because we did not ask for it in the advanced options section.
- Assembled transcripts - This is the one we care most about.



Rename the two **assembled transcripts** datasets to include **pre-pro-B** or **pro-B** in the dataset name. Then add the assembled transcripts datasets to the visualization you created earlier.

Once in Trackster, add the other assembled transcripts dataset and hide the content of the two mapped reads datasets, and zoom out a few times. Set the display mode for both new datasets to Squish.



A couple of things to note about this region:

- the pro-B assembled transcripts have some transcripts that are not in the reference, and even some new exons
- the pre-pro-b assembled transcripts are largely a subset of the pro-B predicted transcripts.
- several annotated transcripts and exons are not detected at all in these datasets.


Save the visualization and go back to **Analyze Data**.

And there is more!

Workflows: repeating an analysis

We'll want to apply the same analysis to our full datasets, and maybe with other similar experiments as well. To do this create a workflow - a repeatable recipe for doing an analysis.

Galaxy supports *de novo* workflow creation, and creating them from histories. To create a workflow from our history, *click* the **cog** at the top of the history panel, and then *select* **Extract Workflow**.



The following list contains each tool that was run to create the datasets in your current history. Please select those that you wish to include in the workflow.

Tools which cannot be run interactively and thus cannot be incorporated into a workflow will be shown in gray.

Workflow name

Tool	History items created
Unknown <i>This tool cannot be used in workflows</i>	1: GTF Annotation chrX/12 regions ✓ Treat as input dataset
Unknown <i>This tool cannot be used in workflows</i>	2: pre-pro-B chrX/12 regions FASTQ raw ✓ Treat as input dataset
Unknown <i>This tool cannot be used in workflows</i>	3: pro-B chrX/12 regions FASTQ raw ✓ Treat as input dataset
FastQC ✓ Include "FastQC" in workflow	4: FastQC on data 2: Webpage
FastQC ✓ Include "FastQC" in workflow	5: FastQC on data 2: RawData
FastQC ✓ Include "FastQC" in workflow	6: FastQC on data 3: Webpage
FastQC ✓ Include "FastQC" in workflow	7: FastQC on data 3: RawData
Trimmomatic ✓ Include "Trimmomatic" in workflow	8: Trimmomatic on pre-pro-B chrX/12 regions FASTQ raw

This takes you to the create workflow page, where you can specify which of the steps and input in your current history you wish to include in the new workflow.

Give the new workflow an informative name and *click* **Create Workflow**. This message appears:



Workflow "RNA-Seq 2 condition QC through Transcript Prediction" created from current history. You can [edit](#) or [run](#) the workflow.

Test the new workflow on the exact same inputs that we just analyzed manually

Click the **run** link in the message (or you can also get there by clicking Workflow in the top bar). This brings up a form asking you to define inputs to the workflow. Set the first three input datasets for this run to the first three datasets from your current history:


Scroll down to the bottom of the form and check **Send results to a new history**. Name the workflow and click **Run workflow**.

This displays an enormous green box in the middle panel. Click the link near the top of the box to get to the new history.

The screenshot shows the Galaxy workflow runner interface for a workflow titled "RNA-Seq 2 condition QC through Transcript Prediction". At the top right are "Expand All" and "Collapse" buttons. The interface is divided into three steps, each with an "Input Dataset" section. Step 1: "Input dataset" shows "1: GTF Annotation chrX/12 regions" selected. Step 2: "Input dataset" shows "2: pre-pro-B chrX/12 regions FASTQ raw" selected. Step 3: "Input dataset" shows "3: pro-B chrX/12 regions FASTQ raw" selected. Each selection is in a dropdown menu with a "type to filter" input field below it.

Working with multiple histories

Running the workflow and sending results to a new history means that we now have two histories in this Galaxy instance. There are a couple of ways to see all your existing histories and to switch between them. The most recent and easiest to use is the all histories view. This can be accessed by clicking on the **table icon (View all histories)** at the top of the history panel.



The all histories view presents all your saved histories (we only have two), with the current one pinned at the left, and your other histories listed in reverse chronological order, from left to right.

Sharing and publishing

Galaxy histories, workflows, and visualization can all be shared and published with Galaxy. Sharing in Galaxy means sharing something with someone else either directly with their Galaxy account, or by creating a URL that can be shared. In addition Galaxy objects can be published, making it easy for anyone to discover the object. Anything that is published will be listed in the appropriate Shared Data section.

Galaxy pages

The datasets and analysis used in the paper (and in this tutorial) are available here:

<https://usegalaxy.org/u/thereddylab/p/prediction-of-gene-activity-in-early-b-cell-development-based-on-an-integrative-multi-omics-analysis>

This is a *Galaxy Page*, a document integrated into a Galaxy server that embeds and provides direct links to Galaxy objects such as histories, workflows, and datasets. Galaxy Pages are a way to bundle together all related analysis and to describe the semantics of your analysis.